# 3M™ Clean-Trace™ ATP荧光检测系统

See Clean  
with Confidence

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4.11.23	As a minimum for food contact surfaces and processing equipment limits of acceptable and unacceptable cleaning performance shall be defined. This shall be based on the potential hazards (e.g. microbiological, allergen, or foreign body contamination or product to product contamination). Acceptable levels of cleaning may be defined by visual appearance, ATP bioluminescence techniques (see Glossary), microbiological testing or chemical testing as appropriate. Where cleaning procedures are part of a defined pre-requisite plan to control the risk of a specific hazard the cleaning and disinfection procedures and frequency shall be validated and records maintained.
4.11.34	The resources for undertaking cleaning shall be available. Where it is necessary to dismantle equipment for cleaning purposes or to enter large equipment for cleaning, this shall be appropriately scheduled and, where necessary, planned for non-production periods. Cleaning staff shall be adequately trained or engineering support provided where access within equipment is required for cleaning.
4.11.45	The cleanliness of equipment shall be checked before equipment is released back into full-production. The results of checks on cleaning, including visual, analytical and microbiological checks, shall be recorded and used to identify trends in cleaning performance and instigate improvements where required.
4.11.66	Cleaning equipment shall be: <ul style="list-style-type: none"> <li>hygienically designed and fit for purpose</li> <li>suitably identified for intended use, e.g. colour coded or labelled</li> <li>cleaned and stored in a hygienic manner to prevent contamination.</li> </ul> Equipment used for cleaning in high-care and high-risk areas shall be visually distinctive and shall be dedicated for use in that area.
4.11.67	Cleaning in place (CIP)
4.11.67.1	Cleaning-in-place (CIP) facilities, where used, shall be monitored and maintained to ensure their effective operation.
4.11.67.2	A schematic plan diagram of the layout of the CIP system shall be available. There shall be an inspection report or other verification-validation that: <ul style="list-style-type: none"> <li>systems are hygienically designed with no dead areas, limited interruptions to flow streams and good system drain ability.</li> <li>scavenge pumps are operated to ensure that there is no build-up of cleaning fluids in the vessels.</li> <li>spray balls effectively clean vessels by providing full surface coverage and are periodically inspected for blockages. Rotating spray devices should have a defined operational time.</li> </ul>

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### Evaluation

### ABSTRACT

practitioners have used the standard method for the detection of hands, the natural ATP biological fluorescence intensity) the value which pick shows that r=0.9

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### Evaluation on ATP bioluminescence assay for rapid determination for Bacterial contamination of water

**ABSTRACT Objective:** To explore the feasibility of the ATP bioluminescence assay for water pollution detection, and providing the scientific basis to water quality monitoring and rapid detection of food safety. **Methods:** At the same time on the Pearl River water and E.coli suspension solution to doing ATP biological fluorescence detection experiment and doing national standard method (Ager plate count method) experiment, the results were compared. **Results:** The higher RLU (relative light unit) value of the ATP biological fluorescent assay obtained, the more total number of live bacteria have, so the ATP biological fluorescent assays with the relevant the national standards method. **Conclusions:** ATP bioluminescence assay can be used on-site rapid detection of water quality monitoring and pollution assessment.

**KEY WORDS** ATP bioluminescence assay, Water quality, The national standard method

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



Extractant  
Left by Swab's  
First Pass Releases  
ATP from Biofilm

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
## What makes a good ATP hygiene monitoring system?

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by Chris Griffith, Von Holy Consulting, Welgemoed, South Africa.

None in the food industry doubts the value of effective cleaning or the role that Adenosine Triphosphate (ATP) bioluminescence can play in monitoring and verifying cleaning protocols. Clean food contact and environmental surfaces are increasingly recognised to be important in limiting pathogen spread, preventing cross contamination and in the manufacture of food and beverage products that are safe and of a quality to meet shelf life needs or claims. The problem is how to decide what makes a good ATP bioluminescence detection system.

### Choosing a system

Marketing and advertising materials should help – but there is a famous quote “there are lies, damn lies and statistics”. One definition of quality

Each ATP test system manufacturer tries to establish unique selling features for their own system. They may highlight individual aspects or advantages of their test or instrument without presenting an overall picture or in relation to industry needs and this can be misleading.

There is relatively little reliable information on what industry actually wants from an ATP hygiene monitoring test system although research reporting on this has subsequently been found to be still valid today (see Table 1).

This published research sets out the role of ATP hygiene monitoring

Attributes/ characteristics	Industry ranking
Speed of results	1
Accuracy of results	2
Reproducibility of results	3
Simplicity and ease of use	4
Cost per test	5
Reliability of equipment	6

important if a large number of tests need to be undertaken or, as is often the case, there is only a short window of opportunity for testing cleaning efficacy before production resumes.

A difference of even 10 seconds per test can add up if several hundred tests need to be performed. To an extent the time taken for a result is going to be dependent upon the time the manufacturer has set for the instrument to read the test and the number of steps or stages required to activate the test.

Collectively this is known as the ‘time to result’ (TTR) which varies by over 60% between some leading instruments. Having a single step in swab activation, apart from being quicker, also helps to reduce variability in results and improve repeatability and reproducibility as well as ease of use.

Ranked 2 and 3 in importance by industry are accuracy and reproducibility and it is perhaps in this area that there is most confusion, claims

consistency of the results between two or more operatives. Clearly both are important, especially when sampling is undertaken by more than one person and results are being used correctly and proactively as part of trend analysis, which is vital to the effective management of cleaning.

Some studies report on factors that actually make no difference to the reliability of the result for real tests on surfaces and are misleading. For example comparing the signal from ATP added to the sampling swab and the signal from ATP that is added to the test reagent which can only be done by dismantling the device, then measuring the difference and reporting this as a measure of ‘precision’ or ‘accuracy’.

In reality, what is important in repeatability and accuracy is the sum of all the components in the test system (see Table 2) and how these factors all interact to achieve the overall final result.

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